

Fisherbrand Disposable Filterware

Increased cell culture sample protection

Integral PES Filter Units and PES Bottle Top Filters

Increase cell culture sample protection with these Fisherbrand™ sterile disposable vacuum PES filter units (with integral receiver bottles) or separate bottle top filters also with PES membranes. Polyethersulfone (PES) membranes have low extractables and low protein binding for use in sterile filtration and clarification of cell culture media and serum.

Integral Filter Unit and Bottle Top Features:

- PES membranes are low protein binding to maintain protein balance
- PES membranes have low extractables to maintain media purity
- Upper chamber of the integral filter units is clearly marked with pore size, membrane type, catalogue number, lot number and expiration date
- Non-pyrogenic and non-cytotoxic
- Gamma radiation-sterilised and individually bagged for a five-year sterile shelf life
- 0.2µm filters for sterile filtration
- 0.45µm filters for fluid clarification and particle removal
- Separate bottle top filters have screw threads for attachment to bottles with 45mm neck sizes



Fisherbrand Storage Receiver Bottles

Reduce the risk of cross contamination with these sterile, disposable polystyrene media bottles. Use in combination with the separate bottle-top filters or as replacement/extra receptacles for the combined filter units. Also excellent as stand-alone single-use storage containers.

Storage Receiver Bottle Features:

- Headspace allows for extra serum or other additives to be added to the filtrate
- Sterilised by gamma radiation
- Individually bagged for a five-year sterile shelf life
- Non-pyrogenic and non-cytotoxic



Cat. No	Description	Capacity, mL	Pore Size (µm)	Pack qty
Fisherbrand Bottle Top Filters				
15953307	Bottle top unit, PES, 50mm filter diameter	150	0.20	12
15963307	Bottle top unit, PES, 50mm filter diameter	150	0.45	12
15973307	Bottle top unit, PES, 75mm filter diameter	500	0.20	12
15983307	Bottle top unit, PES, 75mm filter diameter	500	0.45	12
15993307	Bottle top unit, PES, 90mm filter diameter	1000	0.20	12
Fisherbrand Storage Receiver Bottles				
15903317	Storage receiver bottle, PS, 45mm neck	150	N/A	24
15913317	Storage receiver bottle, PS, 45mm neck	250	N/A	24
15923317	Storage receiver bottle, PS, 45mm neck	500	N/A	12
15933317	Storage receiver bottle, PS, 45mm neck	1000	N/A	12
Fisherbrand Integral Filter Units				
15973297	Filter unit, PES, 50mm filter diameter	150	0.20	12
15983297	Filter unit, PES, 50mm filter diameter	150	0.45	12
15993297	Filter unit, PES, 50mm filter diameter	250	0.20	12
15903307	Filter unit, PES, 50mm filter diameter	250	0.45	12
15913307	Filter unit, PES, 50mm filter diameter	500	0.20	12
15923307	Filter unit, PES, 50mm filter diameter	500	0.45	12
15933307	Filter unit, PES, 90mm filter diameter	1000	0.20	12
15943307	Filter unit, PES, 90mm filter diameter	1000	0.45	12

Find out more at eu.fishersci.com

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
 Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
 Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500



© 2019 Thermo Fisher Scientific Inc. All rights reserved.
 Trademarks used are owned as indicated at fishersci.com/trademarks.